



Technical Data Sheet

CSM2929 CS Easy Prone Standard Covers



Description

- Disposable covers for the CS Easy Prone (all sizes)
- Helps protect the CS Easy Prone from damage or contamination
- Latex and phthalate free
- Single-use aids infection control

Application Area

- OR
- Anaesthesia

Measurements

Chest

- Length: 500mm/19.69"
- Width: 250mm/9.84"
- Height: 170mm/6.69"

Pelvis

- Length: 500mm/19.69"
- Width: 335mm/13.19"
- Height: 170mm/6.69"

Material Composition

- Blue non-woven (waterproof): 55% Woodpulp
45% Polyester
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2929	CS Easy Prone Standard Covers	6 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- If incorrectly applied to the CS Easy Prone, this product can cause pressure marks or ulcers.
- Do not reuse
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624129292

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark