



Technical Data Sheet

CSM2996 CS FlexiFin™ Stirrup Pads (Pair)



Order Information

Product Code	Product Name	Qty
CSM2996	CS FlexiFin™ Stirrup Pads (Pair)	1 Pair

Packaging Information

- Each pair is packaged in a recyclable polyethylene bag
- 1 pair is packaged into a cardboard box for sale

Description

- Reusable - clean using non-alcohol hospital wipes
- Secures and cushions patient's feet and legs in the CS FlexiFin™ Stirrups
- For use in a variety of surgical procedures
- Reduces the incidence of pressure ulcers and marking
- Maximum patient weight: 180kg
- Silicon straps supporting improved infection control

Application Area

- OR

Measurements

- Length mm 800/31.5"
- Width 650mm/25.5"
- Height 20mm/0.8"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- Transfer coated PU with anti-bacterial treatment on knitted polyester
- Silicon

Cleaning

- Clean using non-alcohol hospital wipes

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624129964

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.02 01/04/2025