



# The CS Supine Bariatric Support

A 2-piece wedge cushion set to elevate the head and thorax in the supine position, reducing pressure and facilitating intubation.

- + Helps to easily reach the optimal intubation position by raising the head, upper body and shoulders above the chest ("ramped position")
- + Base can be used alone to allow for optimal positioning for patients of varying size and stature
- + Fully sealed for easy cleaning and improved infection control
- + Fully Radiolucent
- + Maximum Weight: 250kg




## CS Supine Bariatric Support

Code: CSM2648 CS Supine Bariatric Support (Set)  
(Includes Top & Base) Qty: 1

Code: CSM2914 CS Supine Bariatric Support - Top Qty: 1

Code: CSM2915 CS Supine Bariatric Support - Base Qty: 1

## Our recommended addition:



### CS Supine & Lateral Head Support

Code: CSM2643  
Qty: 1 Unit



**Contact us for more information:**  
[www.care-surgical.com](http://www.care-surgical.com) | [info@care-surgical.com](mailto:info@care-surgical.com)