



# Technical Data Sheet

## CSM-2648 CS Supine Bariatric Support



### Order Information

Product Code	Product Name	Qty
CSM-2648	CS Supine Bariatric Support	1 Set

### Description

- 2-part support (CSM-2914 & CSM-2915) to elevate the head and thorax in the supine position, reducing pressure and facilitating intubation.
- Use the base cushion on its own or together the top cushion depending on the stature of the patient
- Reduces the incidence of pressure ulcers and marking
- Reusable

### Application Area

- OR
- Anaesthesia

### Measurements

- Bottom wedge: L 715mm/28.14" | W 515mm/20.27" | H 150mm/5.9"
- Upper wedge: L 785mm/30.9" | W 515mm/20.27" | H 80mm/3.14"

### Material Composition

- Polyurethane foam
- Latex & phthalate free
- Weft knitted polyamide, coated with anti-microbial polyurethane coating

### Cleaning

- Clean using alcohol-based hospital wipes

### Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 1 set is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Do not store anything heavy on top
- Avoid contact with sharp objects

### Contraindications

- Incorrect positioning of the patient on the CS Supine Bariatric Support can cause:
- Neurological damage, which can cause paraplegia;
- Venous pressure, which can lead to ischaemia.
- Pressure marking or ulcers
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel
- Do not immerse in water
- Do not autoclave

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No:

Top (CSM-2914) - 5060624129148

Base (CSM-2915) - 5060624129155

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.02 01/03/2024