



Technical Data Sheet

CSM2699 CS Easy Prone™ - 150mm - Chest & Pelvis Support



Description

- Set consists of: CSM2920 CS Easy Prone™ - 150mm - Chest Only & CS2921 CS Easy Prone™ - 150mm - Pelvis Only
- Raises the patient's chest and pelvis in the prone position
- Enables compression-free positioning of the abdomen
- Reusable, fully radiolucent and fully sealed
- The design prevents pressure in the armpit, neck and genital area
- Suitable for many different body shapes, as it can be flexibly adjusted
- Enables a natural lordosis of the spine
- Allows movements of the table without compromises
- Weight capacity of 150kg

Application Area

- OR
- Anaesthesia

Measurements

Chest:

- Length 500mm/19.69"
- Width 250mm/9.84"
- Height 150mm/5.91"

Pelvis:

- Length 500mm/19.69"
- Width 336mm/13.23"
- Height 150mm/5.91"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- Weft knitted polyamide, coated with anti-microbial polyurethane coating

Order Information

Product Code	Product Name	Qty
CSM2699	CS Easy Prone™ - 150mm - Chest & Pelvis Support	1 Set

Packaging Information

- Each component of the set is packaged in a recyclable polyethylene bag
- 1 set is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place

Contraindications

- Incorrect positioning of the patient on the CS Easy Prone™ can cause:
 - Pressure marking or ulcers
 - Neurological damage, which can cause paraplegia;
 - Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624126994

Cleaning

- Clean using alcohol-based hospital wipes

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device Reuse Latex Free Non Sterile UK Conformity Assessed CE Mark

Rev.03 05/03/2026