



Technical Data Sheet

CSM2918 CS Easy Prone™ - 170mm - Chest & Pelvis Support



Description

- Set consisting of CSM2922 CS Easy Prone™ - 170mm - Chest Only & CSM2923 CS Easy Prone™ - 170mm - Pelvis Only
- Raises the patient's chest and pelvis in the prone position
- Enables compression-free positioning of the abdomen
- Reusable, fully radiolucent and fully sealed
- The design prevents pressure in the armpit, neck and genital area
- Suitable for many different body shapes, as it can be flexibly adjusted
- Enables a natural lordosis of the spine
- Allows movements of the table without compromises
- Weight capacity of 160kg

Application Area

- OR
- Anaesthesia

Measurements

Chest:

- Length 500mm/19.69"
- Width 270mm/10.63"
- Height 170mm/6.69"

Pelvis:

- Length 500mm/19.69"
- Width 356mm/14.02"
- Height 170mm/6.69"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- Weft knitted polyamide, coated with anti-microbial polyurethane coating

Order Information

Product Code	Product Name	Qty
CSM2918	CS Easy Prone™ - 170mm - Chest & Pelvis Support	1 Set

Packaging Information

- Each component of the set is packaged in a recyclable polyethylene bag
- 1 set is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place

Contraindications

- Incorrect positioning of the patient on the CS Easy Prone™ can cause:
- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia;
- Venous pressure, which can lead to ischaemia.
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624129186

Cleaning

- Clean using alcohol-based hospital wipes

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
M43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Rev.04 05/03/2026