



Technical Data Sheet

CSM2919 CS Easy Prone™ - 200mm - Chest & Pelvis Support



Description

- Set consists of: CSM2924 CS Easy Prone™ - 200mm - Chest Only & CSM2925 CS Easy Prone™ - 200mm - Pelvis Only
- Raises the patient's chest and pelvis in the prone position
- Enables compression-free positioning of the abdomen
- Reusable, fully radiolucent and fully sealed
- The design prevents pressure in the armpit, neck and genital area.
- Suitable for many different body shapes, as it can be flexibly adjusted
- Enables a natural lordosis of the spine
- Allows movements of the table without compromises
- Weight capacity of 180kg

Application Area

- OR
- Anaesthesia

Measurements

Chest:

- Length: 500mm/19.69"
- Width: 290mm/11.42"
- Height: 200mm/7.87"

Pelvis:

- Length: 500mm/19.69"
- Width: 376mm/14.8"
- Height: 200mm/7.87"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- Weft knitted polyamide, coated with anti-microbial polyurethane coating

Order Information

Product Code	Product Name	Qty
CSM2919	CS Easy Prone™ - 200mm - Chest & Pelvis Support	1 Set

Packaging Information

- Each component of the set is packaged in a recyclable polyethylene bag
- 1 set is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place

Contraindications

- Incorrect positioning of the patient on the CS Easy Prone™ can cause:
 - Pressure marking or ulcers
 - Neurological damage, which can cause paraplegia;
 - Venous pressure, which can lead to ischaemia.
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624129193

Cleaning

- Clean using alcohol-based hospital wipes

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Ref: CSM2919-TDS-(Rev.04) 05/03/2026