



CARE SURGICAL


CS Supine Bariatric Support

Product Code: CSM2648



Only to be used by trained clinicians. Ensure product is secure before positioning patient.

The CS Supine Bariatric Support reduces pressure on the chest and facilitates intubation.



2-piece support set for elevated positioning of the patient's head and thorax in the supine position.

Can be used with Top (CSM2914) and Base (CSM2915) or only with Base.


Cleaning instructions:



Clean using standard hospital wipes



Do not immerse in water



Do not autoclave



UK: Care Surgical Ltd
Unit 6 Ringtail Road,
Burscough,
L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



SAFE DISPOSAL: customers should adhere to any local laws and regulations to the safe disposal of medical devices.

DO NOT USE products that appear to be damaged or worn. In an event of product damage or patient injury, please contact Care Surgical.



5 060624 120381



Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland

Ref: CSM2648 - IFU - (Rev.04) 04/2025