



Technical Data Sheet

CSM2754 CS FlexiFin™ Stirrups (Pair)



Description

- Light stirrup leg holders to protect the patients calves in lithotomy position.
- Includes silicone strap pad attachment system - CS FlexiFin™ Stirrup Pads
- Reusable - clean using alcohol-based hospital wipes
- Latex and phthalate free
- Patient Weight Capacity: 227kg/500lbs.

Application Area

- OR

Measurements

- Length 1000mm/39.37"
- Width 680mm/26.77"
- Height 560mm/22"

Material Composition

Stirrup:

- Stainless Steel
- Aluminium
- ABS (Boot)
- Rubber (Shrouds and housings)
- Hook velcro

Pads:

- Polyurethane coated fabric
- Viscoelastic memory foam
- Silicone
- Loop velcro

Cleaning

- Clean using non-alcohol hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2754	CS FlexiFin™ Stirrups (Pair)	1 Pair

Packaging Information

- Each pair is packaged in a recyclable polyethene bag and bubblewrapped
- 1 pair is packaged into a cardboard box for sale

Storage Instructions

- Transport and store in dry conditions

Contraindications

- Incorrect attachment to the surgical tables could result injuries to the patient or nurse
- Do not use if the blade clamp appears to be defective or damaged
- Do not use if the Stirrup Legholders appear to be defective or damaged
- Do not autoclave
- Only to be used by trained medical personnel
- Do not exceed the device safe working load
- Do not immerse in water
- After each use, clean the device with alcohol-based wipes
- Make sure that the device is dry before you store it or use it again

Features

- Patient Weight Capacity: 227kg/ 500lbs.
- Lithotomy adjustment (from sterile field): Yes
- Lithotomy range (from sterile field): +90° to -45°
- Abduction range (from sterile field): +25° to -10°
- Maximum medial/lateral shift/ boot rotational angle: +-23°
- Product weight: Approx. 13kg per pair

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624127540
- Warranty:
 - 2-year warranty on stirrup components
 - 1-year warranty on boot pads

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor:

MDSS AU Pty Ltd, Ground Floor, Suite 55,
97-99 Bathurst Street, Sydney, NSW 2000, Australia

Rev.07 05/03/2026