



## The CS Carbon Spine Extension - Merivaara

The CS Carbon Spine Extension can be fitted to a range of standard operating tables offering a patient position and imaging capabilities comparable to that of a specialist dual-rail spine table.



- + Improved efficiency- simply remove the extension to return the table for standard use
- + Greatly improved anaesthetic access- no end columns between the anaesthetic team and the patient
- + 360 degree C-Arm/0-Arm access
- + Significantly more cost-effective in comparison to purchasing a full spine table
- + Detachable support leg for added stability if preferred
- + Light weight – one person required to set up
- + 250kg weight capacity

### CS Carbon Spine Extension - Merivaara

Code: CSM-2930  
Qty: 1 Set

Contains:  
CSM-2730 - CS Carbon Spine Extension Frame

CSM-2737 - CS Carbon Spine Extension Leg

CSM-2742 - CS Carbon Spine Extension Adaptor - Merivaara



### Related Products



**CS Prone Plus Head Support System**

Code: CSM-2565  
Qty: 1 System  
(Helmet & Mirror)



**CS Universal Arm support - Moulded**

Code: CSM-2647  
Qty: 1 Unit



**CS Shin Supports (Pair)**

Code: CSM-2820  
Qty: 1 Pair



**CS Knee Rings (Pair)**

Code: CSM-2821  
Qty: 1 Pair

Find Out More:



**Contact us for more information:**  
[www.care-surgical.com](http://www.care-surgical.com) | [info@care-surgical.com](mailto:info@care-surgical.com)