



# Technical Data Sheet

## CSM-2730 CS Carbon Spine Extension Frame



### Description

- A dual rail extension intended to be attached to an operating table, to which specialist supports can be attached to position patients for surgery (e.g. spinal surgery)
- Compatible with standard operating tables
- Compatible with sectional pressure management pads and other support attachments (hip/thigh supports)
- Reusable

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 1250mm/49.21"
- Width 540mm/21.26"
- Height 64mm/2.52"

### Material Composition

- Carbon Fibre
- Stainless Steel 316
- Aluminium 6061
- Cerakote
- Loctite
- 3M Scotchweld EPX DP-490
- Synthetic Rubber TBD

### Cleaning

- Clean using non-alcohol hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM-2730	CS Carbon Spine Extension Frame	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

### Contraindications

- Ensure the platform is securely attached to the table before use
- Check the platform for damage before and after each use
- Do not use if damaged Contact supplier for advice
- Do not exceed a maximum weight of 250kg
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624127304
- Latex and phthalate free

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson,  
MA 01749, USA  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 02/04/2024



# Technical Data Sheet

## CSM-2737 CS Carbon Spine Extension Leg



### Description

- Reusable - clean using non-alcohol hospital wipes.
- Only for use with the CS Carbon Spine Extension
- Adds stability to the CS Carbon Spine Extension
- Helps to reduce movement of the CS Carbon Spine Extension

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 80mm/3.15"
- Width 120mm/4.72"
- Height 1200mm/47.24"

### Material Composition

- Stainless Steel (304 & 316)
- Aluminium 6061
- Nylon 66
- Onyx Filament
- ABS Plastic
- PVC Rubber
- EPDM Rubber

### Cleaning

- Clean using non-alcohol hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM-2737	CS Carbon Spine Extension Leg	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

- Incorrect attachment to the CS Carbon Spine Extension could result in movement of the spine extension frame.
- Do not use if the leg appears to be defective or damaged.
- Only to be used by trained medical personnel
- Maximum weight limit: 250kg
- Do not immerse in water
- Do not autoclave

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624123320
- Latex and phthalate free

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson,  
MA 01749, USA  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 02/04/2024



# Technical Data Sheet

## CSM-2742 CS Carbon Spine Extension Adaptor - Merivaara



### Description

- Adaptor to connect the CS Carbon Spine Extension frame to Merivaara operating tables
- Compatible with Merivaara operating table
- Fixed position enabled by static table fixtures
- Reusable - clean using non-alcohol hospital wipes

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 180mm/7"
- Width 540mm/21"
- Height 75mm/3"

### Material Composition

- Stainless Steel
- Aluminium
- Carbon Fibre
- ABS/Nylon Plastic

### Cleaning

- Clean using non-alcohol hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM-2742	CS Carbon Spine Extension Adaptor - Merivaara	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag and bubblewrap, if applicable
- Each unit is packaged into a cardboard box for sale.

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product.

### Contraindications

- Ensure the adaptor is securely attached to the table before use
- Check the adaptor for damage before and after each use. Do not use if damaged
- Contact supplier for advice
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624125096
- Latex and phthalate free

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson,  
MA 01749, USA  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 02/04/2024