



Technical Data Sheet

CSM2730 CS Carbon Spine Extension Frame



Description

- A dual rail extension intended to be attached to an operating table, to which specialist supports can be attached to position patients for surgery (e.g. spinal surgery)
- Compatible with standard operating tables
- Compatible with sectional pressure management pads and other support attachments (hip/thigh supports)
- Reusable

Application Area

- OR
- Anaesthesia

Measurements

- Length: 1250mm/49.21"
- Width: 540mm/21.26"
- Height: 64mm/2.52"

Material Composition

- Carbon Fibre
- Stainless Steel 316
- Aluminium 6061
- Cerakote
- Loctite
- 3M Scotchweld EPX DP-490
- Synthetic Rubber TBD
- Latex and phthalate free

Order Information

Product Code	Product Name	Qty
CSM2730	CS Carbon Spine Extension Frame	1 Unit

Cleaning

- Clean using non-alcohol hospital wipes

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

- Ensure the platform is securely attached to the table before use
- Check the platform for damage before and after each use
- Do not use if damaged Contact supplier for advice
- Do not exceed a maximum weight of 250kg
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624127304

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
443 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark



Technical Data Sheet

CSM2737 CS Carbon Spine Extension Leg



Description

- Reusable - clean using non-alcohol hospital wipes.
- Only for use with the CS Carbon Spine Extension
- Adds stability to the CS Carbon Spine Extension
- Helps to reduce movement of the CS Carbon Spine Extension

Application Area

- OR
- Anaesthesia

Measurements

- Length: 80mm/3.15"
- Width: 120mm/4.72"
- Height: 1200mm/47.24"

Material Composition

- Stainless Steel (304 & 316)
- Aluminium 6061
- Nylon 66
- Onyx Filament
- ABS Plastic
- PVC Rubber
- EPDM Rubber
- Latex and phthalate free

Cleaning

- Clean using non-alcohol hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2737	CS Carbon Spine Extension Leg	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- Incorrect attachment to the CS Carbon Spine Extension could result in movement of the spine extension frame.
- Do not use if the leg appears to be defective or damaged.
- Only to be used by trained medical personnel
- Maximum weight limit: 250kg
- Do not immerse in water
- Do not autoclave

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624123320

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark



Technical Data Sheet

CSM2742 CS Carbon Spine Extension Adaptor - Merivaara



Description

- Adaptor to connect the CS Carbon Spine Extension frame to Merivaara operating tables
- Compatible with Merivaara operating table
- Fixed position enabled by static table fixtures
- Reusable - clean using non-alcohol hospital wipes

Application Area

- OR
- Anaesthesia

Measurements

- Length: 180mm/7"
- Width: 540mm/21"
- Height: 75mm/3"

Material Composition

- Stainless Steel
- Aluminium
- Carbon Fibre
- ABS/Nylon Plastic
- Latex and phthalate free

Cleaning

- Clean using non-alcohol hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2742	CS Carbon Spine Extension Adaptor - Merivaara	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag and bubblewrap, if applicable
- Each unit is packaged into a cardboard box for sale.

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product.

Contraindications

- Ensure the adaptor is securely attached to the table before use
- Check the adaptor for damage before and after each use. Do not use if damaged
- Contact supplier for advice
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624125096

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark



Technical Data Sheet

CSM-2751 CS Head Tray Attachment - Lowered



Order Information

Product Code	Product Name	Qty
CSM2751	CS Head Tray Attachment - Lowered	1 Unit

Description

- To be used with a head support and/or face cushion to support a patient's head
- Secure the tray to the side bars of the operating table and lock onto place with the quick-release locking clips
- Reusable - clean using non-alcohol hospital wipes

Application Area

- OR
- Anaesthesia

Measurements

- Length: 501mm/20"
- Width: 296mm/12"
- Height: 155mm/6"

Material Composition

- Aluminium
- Stainless Steel
- Carbon Fibre plates
- Latex and phthalate free

Cleaning

- Clean using non-alcohol hospital wipes

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- Ensure the platform is securely attached to the table before use
- Check the platform for damage before and after each use
- Do not use if damaged Contact supplier for advice
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624127519

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark



Technical Data Sheet

CSM2560 CS Prone Head Support System



Description

- Reusable
- Helmet holds the CS Prone Face Cushion (CSM2510) to maintain their integrity throughout prone procedures
- Space for safe positioning of ET tubes
- The helmet can be used with or without the mirror
- The mirror enables the clinician to monitor the patient's eyes and tubing during the procedure
- The cogs fitted to the mirror are used to adjust the height or angle of the helmet to achieve the safest natural neutral neckline

Application Area

- OR
- Anaesthesia
- ICU

Measurements

- Helmet:
 - Length: 310mm/12.2"
 - Width: 220mm/8.6"
 - Height: 90mm/3.5"
- Mirror:
 - Length: 220mm/8.6"
 - Width: 220mm/8.6"
 - Height: 4mm/0.16"

Material Composition

- Helmet made from ABS Plastic
- Acrylic mirror, which incorporates 4 adjustable cogs
- Aluminium 6061
- Latex and phthalate free

Order Information

Product Code	Product Name	Qty
CSM2560	CS Prone Head Support System	1 System

Cleaning

- Clean with standard hospital wipes

Packaging Information

- Each set is packaged in a recyclable polyethylene bag
- Each set is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- POVL
- pressure ulcers
- neurological damage
- venous pressure, which can lead to ischaemia
- Do not use if there is any damage on the helmet or mirror
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Not fully radiolucent- cog screws may be visible on x-rays
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624120114

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark



Technical Data Sheet

CSM2510 CS Prone Face Cushion



Description

- Prone Face Cushion - Made from polyurethane foam
- For use with the CS Head Support System and other compatible head support systems
- Distributes pressure on the patient's cheeks and forehead while in the prone position
- For patients aged 8 and over, dependent on stature
- For patients weighing up to 200kg

Application Area

- OR
- Anaesthesia

Measurements

- Length: 312mm/12.28"
- Width: 224mm/8.82"
- Height: 146mm/5.75"

Material Composition

- Polyurethane foam
- Latex and phthalate free

Order Information

Product Code	Product Name	Qty
CSM2510	CS Prone Face Cushion	10 Units

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 10 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- POVL
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624120015

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not
Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark



Technical Data Sheet

CSM2647 CS Universal Arm Support - Moulded



Description

- Supports the patient's arm when undergoing surgical procedures
- Reusable
- Attaches to most standard and specialist operating tables to support patients' arms in various positions during surgical procedures
- Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)
- Supports maximum weight of 10kg

Application Area

- OR
- Anaesthesia

Measurements

- Length: 460mm/18.11"
- Width: 150mm/5.91"
- Height: 550mm/21.65"

Material Composition

- Stainless Steel and PU Foam
- Latex and phthalate free

Cleaning

- Clean using non-alcohol hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2647	CS Universal Arm Support - Moulded	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624123511

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark



Technical Data Sheet

CSM2739-1 CS Table Rail Adaptor EU



Description

- Adaptor for use on standard dual-rail tables, with EU table rail. It is used to provide a secure method of mounting (holding) other medical devices, such as armboards
- Reusable - clean using alcohol-based hospital wipes
- For use with dual rail tables, compatible with CS Carbon Spine Extension Frame

Application Area

- OR
- Anaesthesia

Measurements

- Length 155mm/6.10"
- Width 85mm/3.35"
- Height 70mm/2.76"

Material Composition

- Aluminium,
- Stainless steel
- Plastic
- Latex and phthalate free

Cleaning

- Clean using alcohol-based hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2739-1	CS Table Rail Adaptor EU	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place

Contraindications

- Incorrect attachment to dual rail tables or the CS Carbon Spine Extension Frame could result in the attachments failing
- Do not use if the adaptor appears to be defective or damaged
- Only to be used by trained medical personnel
- Do not exceed the device safe working load
- Do not immerse in water
- Do not autoclave
- After each use, clean the device with alcohol-based wipes
- Make sure that the device is dry before you store it or use it again

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624123573

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark



Technical Data Sheet

CSM2709 CS Chest Support



Description

- Supports the patient's chest in the prone position
- Attach to the CS Spine Frame or other suitable modular operating frame or table to support the patient's chest whilst in the prone position
- Use the incorporated clips to attach the support to the table
- Open the clip to slide the chest support laterally into the required position
- Close the clip securely to ensure the support does not move during the procedure

Application Area

- OR
- Anaesthesia

Measurements

- Length: 208mm/8.19"
- Width: 480mm/18.90"
- Height: 175mm/6.89"

Material Composition

- Wipeable skin covering
- Polyurethane foam
- Transfer coated PU with anti-bacterial treatment on knitted polyester
- Latex and phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2709	CS Chest Support	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel
- If the clips are not fully closed, the support could move during the procedure which could cause injury to the patient or clinician

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624120046

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark



Technical Data Sheet

CSM2717 CS Hip Support (Pair)



Description

- Supports the patient's hip in the prone position
- Attach to the CS Spine Frame or other suitable modular operating frame or table to support the patient's hip(s) whilst in the prone position
- Use the incorporated clip to attach the support to the table
- Open the clip to slide the chest support laterally into the required position
- Close the clip securely to ensure the support does not move during the procedure

Application Area

- OR
- Anaesthesia

Measurements

- Length: 230mm/9.06"
- Width: 136mm/5.35"
- Height: 242mm/9.53"

Material Composition

- Wipeable skin covering
- Polyurethane foam
- Transfer coated PU with anti-bacterial treatment on knitted polyester
- Latex and phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2717	CS Hip Support (Pair)	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel
- If the clips are not fully closed, the support could move during the procedure which could cause injury to the patient or clinician

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624120756

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark



Technical Data Sheet

CSM2821 CS Knee Rings (pair)



Description

- Supports the patients' knees in the prone position
- Reduces the risk of pressure ulcers

Application Area

- OR
- Anaesthesia

Measurements

- Length: 230mm/9.05"
- Width 230mm/9.05"
- Height: 60mm/2.36"

Material Composition

- Polyurethane foam core
- Transfer coated PU with anti-bacterial treatment on knitted polyester
- Latex and phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2821	CS Knee Rings	1 Pair

Packaging Information

- Each pair is packaged in a recyclable polyethylene bag
- 1 pair is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia;
- Venous pressure, which can lead to ischaemia.
- Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624121128

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed




CE Mark



Technical Data Sheet

CSM3020 CS Easy Shin Support™



Description

- Reusable
- Supports the patient's shins in the prone position and allows their ankles/feet to find their own position
- Reduces the incidence of pressure-related injury to the shins

Application Area

- OR
- Anaesthesia

Measurements

- Length: 500mm/19.69"
- Width 220mm/8.66"
- Height: 240mm/9.45"

Material Composition

- Polyurethane foam core
- PU Cover
- Latex & phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM3020	CS Easy Shin Support™	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can lead to temporary or permanent paralysis
- Venous pressure, which can lead to ischaemia
- Examine the product for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624125966

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor:

MDSS AU Pty Ltd, Ground Floor, Suite 55,
97-99 Bathurst Street, Sydney, NSW 2000, Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark