



Technical Data Sheet

CSM3084 CS Carbon Spine Extension Cart



Order Information

Product Code	Product Name	Qty
CSM3084	CS Carbon Spine Extension Cart	1 Unit

Description

- CS Carbon Spine Extension cart made of plastic, stainless steel and aluminium on 4 antistatic, brake-able castors
- Cart that holds, stores and transports 1x CSM2730, CSM2737 & table adaptor along with capacity for a full body positioning set including supports for head, arms, chest, hip, thighs, knees and feet
- Modular design to allow for optimum storage
- Carbon Extension frame and adaptor can be stored on hooks whilst additional supports can be stored on multiple aluminum bars, popper plate and base plate.
- Four antistatic double castors allow easy transport
- All castors can be braked for optimal safety
- Easy to clean with wipe disinfectant

Application Area

- OR

Measurements

- Length 510mm/20.08"
- Width: 600mm/23.62"
- Overall Height: 1720mm/67.7"
- Castors: 125mm/4.92"

Material Composition

- Stainless Steel
- Alu Bars
- HPL
- Rubber

Packaging Information

- Each cart is wrapped with protective foam and packaged in a recyclable polyethene bag
- Each cart is packaged into a cardboard box for sale
- Allen key and screw caps to be sent with item

Storage Instructions

- Transport and store in dry conditions
- Do not store anything heavy on top

Contraindications

- Always hook or clamp supports properly before moving cart to prevent damage.
- Always brake the castors when the cart is not in motion to avoid rolling away and accidents
- Do not damage surfaces with sharp objects
- Do not step on it

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:N
- Basic UDI/GTIN No: 5060624127861

Cleaning

- Clean using non-alcohol hospital wipes

Labelling



Manufacturer




Maximum Load Capacity

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Rev.01 01/07/2026