



Technical Data Sheet

CSM3043 CS Fixed Blade Clamp



Description

- For use on compatible surgical tables in surgical environments and ICUs
- Reusable - clean using alcohol-based hospital wipes
- Latex and phthalate free
- Intended to fix or mount blade-style instrumentation to operating table rails to the device
- Suitable for 25*10mm (EU), 28.5*9.5mm (US), and 30*8mm (JPN) side rails

Application Area

- OR

Measurements

- Length 75mm/2.95"
- Width 45mm/1.77"
- Height 50mm/1.97"

Material Composition

- Latex & phthalate free
- Main Body: - Stainless Steel
- Screw - Stainless Steel & Plastic (ABS)

Order Information

Product Code	Product Name	Qty
CSM3043	CS Fixed Blade Clamp	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Transport and store in dry conditions

Contraindications

- Incorrect attachment to the surgical tables could result in the attachments failing
- Do not use if the blade clamp appears to be defective or damaged
- Only to be used by trained medical personnel
- Do not exceed the device safe working load
- Do not immerse in water
- Do not autoclave
- After each use, clean the device with alcohol-based wipes
- Make sure that the device is dry before you store it or use it again

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624126697

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.02 22/07/2025