






The CS Universal Spine Extension™

The CS Universal Spine Extension™ can be fitted to a wide range of standard operating tables offering a patient position and imaging capabilities comparable to that of a specialist dual-rail spine table.

- + **Lightweight and fully modular** - single person set-up/breakdown
- + **Low profile end column** - unrestricted anaesthetic access
- + **Fully electric height adjustment** - no handle winding
- + **Compatible with Mayfield/Doro Skull Clamps**
- + **Hinge mechanism** to allow flexion and extension of the spine
- + **Small Storage Footprint** - entire system and cart takes up less than 1m space
- + **275kg patient weight capacity**



FIND OUT MORE:



CS Universal Spine Extension™ - EU

Code: CSM2958
Qty: 1 Set



What's Included:

CS Universal Spine Extension™ - Set (EU)

Code: CSM2958
Qty: 1 Set

Includes:
CSM2957 - CS Universal Spine Extension™ - Frame
CSM2955 - CS Universal Spine Extension™ - Table Adaptor (EU)
CSM2956 - CS Universal Spine Extension™ - Powered Leg



CS Prone Plus Head Support System™

Code: CSM2565
Qty: 1 System
(Helmet & Mirror)



CS Universal Arm Support - Moulded™

Code: CSM2647
Qty: 1 Unit



CS Easy Shin Support™

Code: CSM3020
Qty: 1 Unit



CS Knee Rings (Pair)

Code: CSM2821
Qty: 1 Pair



CS Table Rail Adaptor - EU

Code: CSM2739-1
Qty: 1 Unit



CS Prone Plus Face Cushion

Code: CSM2512
Qty: 10 Units



CS Head Tray Attachment - Lowered

Code: CSM2751
Qty: 1 Unit



CS Chest Support

Code: CSM2709
Qty: 1 Unit



CS Hip Supports (pair)

Code: CSM2717
Qty: 1 Pair



CS Thigh Supports (pair)

Code: CSM2718
Qty: 1 Pair



CS Pad 2"

Code: CSM2568-2
Qty: 1 Unit



CS Comfort Covers - 5 Piece


Code: CSM2600
Qty: 6 Units

Optional Accessories:



CS Din Rail Clamp

Code: CSM2763
Qty: 1 Unit



CS Supine Top-Base Only

Code: CSM2722-7
Qty: 1 Unit



CS Supine Top - 2" Foam Pad

Code: CSM2722-2
Qty: 1 Unit



CS Comfort Covers - 5 Piece

Code: CSM2600
Qty: 6 Units



CS Universal Spine Extension™ - Powered Leg Cover

Code: CSM3029
Qty: 25 Units



CS Universal Spine Extension Cart

Code: CSM2954
Qty: 1 Unit

Contact us for more information:
www.care-surgical.com | info@care-surgical.com