






# The CS Universal Spine Extension™

The CS Universal Spine Extension™ can be fitted to a wide range of standard operating tables offering a patient position and imaging capabilities comparable to that of a specialist dual-rail spine table.

- + **Lightweight and fully modular** - single person set-up/breakdown
- + **Low profile end column** - unrestricted anaesthetic access
- + **Fully electric height adjustment** - no handle winding
- + **Compatible with Mayfield/Doro Skull Clamps**
- + **Hinge mechanism** to allow flexion and extension of the spine
- + **Small Storage Footprint** - entire system and cart takes up less than 1m space
- + **600lbs. patient weight capacity**



FIND OUT MORE:



**CS Universal Spine Extension™ - Set (US)**

Code: CSM3031  
Qty: 1 Set



## Whats Included:



### CS Universal Spine Extension™ (Set) - US

Code: CSM3031  
Qty: 1 Set

Includes:

- CSM 2957 - CS Universal Spine Extension™ - Frame
- CSM 3030 - CS Universal Spine Extension™ - Table Adaptor (US)
- CSM 2956 - CS Universal Spine Extension™ - Powered Leg



**CS Prone Plus Head Support System™**

Code: CSM2565  
Qty: 1 System  
(Helmet & Mirror)



**CS Universal Arm Support - Moulded™**

Code: CSM2647  
Qty: 1 Unit



**CS Easy Shin Support™**

Code: CSM3020  
Qty: 1 Unit



**CS Knee Rings (Pair)**

Code: CSM2821  
Qty: 1 Pair



**CS Table Rail Adaptor - US**

Code: CSM2739-3  
Qty: 1 Unit



**CS Prone Plus Face Cushion**

Code: CSM2512  
Qty: 10 Units



**CS Head Tray Attachment - Lowered**

Code: CSM2751  
Qty: 1 Unit



**CS Chest Support**

Code: CSM2709  
Qty: 1 Unit



**CS Hip Supports (pair)**

Code: CSM2717  
Qty: 1 Pair



**CS Thigh Supports (pair)**

Code: CSM2718  
Qty: 1 Pair



**CS Pad 2"**

Code: CSM2568-2  
Qty: 1 Unit



**CS Comfort Covers - 5 Piece**

Code: CSM2600  
Qty: 6 Units

## Optional Accessories:



**CS Din Rail Clamp - US**

Code: CSM2765  
Qty: 1 Unit



**CS Supine Top-Base Only**

Code: CSM2722-7  
Qty: 1 Unit



**CS Supine Top - 2" Foam Pad**

Code: CSM2722-2  
Qty: 1 Unit



**CS Comfort Covers - 5 Piece**

Code: CSM2600  
Qty: 6 Units



**CS Universal Spine Extension™ - Powered Leg Cover**

Code: CSM3029  
Qty: 25 Units



**CS Universal Spine Extension Cart**

Code: CSM2954  
Qty: 1 Unit

Contact us for more information:  
[www.care-surgical.com](http://www.care-surgical.com) | [info@care-surgical.com](mailto:info@care-surgical.com)