




CARE SURGICAL

CS Lateral Support Comfort Covers

Product Code: CSM-3040



Suitable for 3-point fixation in the lateral position (2x iliac crest and 1x sacrum).

Ensures additional comfort through extra padding, enhances infection control and protects the pads from damage. Can be used on both, reusable or disposable pads.

Caution: 100% protection against infection is not guaranteed. Therefore, disposable pads should not be reused.



UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



LOT

EC REP

Ref: CSM-3040 - IFU Rev 01 01/2025