



Technical Data Sheet

CSM3053 CS FlexiFinless™ Stirrup Pads (Pair)



Description

- Reusable - clean using non-alcohol hospital wipes
- Stirrup pads to be used in conjunction with the CS FlexiFinless™ stirrups (Pair) (CSM3050) to protect and support the patients lower limb during surgery
- Reduces the incidence of pressure ulcers and marking
- Silicone straps supporting improved infection control

Application Area

- OR

Measurements

- Length: 630mm/24.8"
- Width: 650mm/25.5"
- Height: 20mm/0.8"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- Transfer coated PU with anti-bacterial treatment on knitted polyester
- Silicone

Cleaning

- Clean using non-alcohol hospital wipes

Order Information

Product Code	Product Name	Qty
CSM3053	CS FlexiFinless™ Stirrup Pads (Pair)	1 Pair

Packaging Information

- Each pair is packaged in a recyclable polyethylene bag
- 1 pair is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624126963

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark