




# CARE SURGICAL

Product Code: CSM3029

## CS Universal Spine Extension- Powered Leg Cover



Transparent cover for CSM2956 CS Universal Spine Extension - Powered Leg to protect the device against liquids.

Place the cover completely over the CSM2956 CS Universal Spine Extension - Powered Leg. If the protruding corners at the top are in the way, they can be folded over and secured with the adhesive dots located at the corners.



UK: Care Surgical Ltd  
Unit 6 Ringtail Road,  
Burscough,  
L40 8JY  
Tel: 01704 336671

USA: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson,  
MA 01749, USA  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,  
Lucan Co. Dublin, K78 X5P8 Ireland



**DO NOT USE** products that appear to be damaged or worn.  
In an event of product damage, please contact Care Surgical.



EU REP