



CARE SURGICAL
Specialists in Patient Positioning

The New CS FlexiFin™ & CS FlexiFin™ Max Lithotomy Stirrups



CS FlexiFin™ Max Stirrups
(Pair) - 360kg/790lbs.

Code: CSM2913

Qty: 1 Pair

CS FlexiFin™ Stirrups
(Pair) - 227kg/500lbs.

Code: CSM2754

Qty: 1 Pair

The CS FlexiFin Stirrups are designed for effortless, optimal lithotomy positioning.



Flexible Boot to accommodate larger patient calves and pressure-relieving padding to reduce the risk of pressure injury



CS FlexiFin™ Stirrups

Other handles

Longer handle for ease of adjustment



Velcro free silicone straps for improved infection prevention



What's included - CS FlexiFin Stirrups CSM2574:



CS FlexiFin™ Stirrups (Pair)

Code: CSM2574
Qty: 1 Pair



CS FlexiFin™ Stirrup Pads (Pair)

Code: CSM2996
Qty: 1 Pair

What's included - CS FlexiFin Stirrups Max CSM2913:



CS FlexiFin™ Max Stirrups (Pair)

Code: CSM2913
Qty: 1 Pair



CS FlexiFin™ Max Stirrup Pads (Pair)

Code: CSM2994
Qty: 1 Pair

Optional additions:



CS Stirrup Cart

Code: CSM2724
Qty: 1 Unit



CS Blade Clamp - UK/EU

Code: CSM2761
Qty: 1 Unit



CS Blade Clamp Plus - EU

Code: CSM2767
Qty: 1 Unit



CS Blade Clamp Plus - US

Code: CSM2762
Qty: 1 Unit

CS FlexiFin & FlexiFin Max Stirrups - Hook & Loop Straps



CS FlexiFin Max Stirrups - Hook & Loop (Pair)

Code: CSM2756
Qty: 1 Pair

CS FlexiFin Stirrups - Hook & Loop (Pair)

Code: CSM2755
Qty: 1 Pair

Features

Patient Weight Capacity	Standard: 227kg/500lbs. Max: 360kg/790lbs.
Lithotomy Adjustments (From sterile field)	Yes
Lithotomy Range (From sterile field)	+90° to -45°
Abduction Range (From sterile field)	+25° to -10°
Maximum medial/lateral shift/boot rotational angle	+/-23°