Technical Data Sheet

CSM3040 CS Lateral Support Comfort Covers



Order Information

Product Code	Product Name	Qty
CSM3040	CS Lateral Support Comfort Covers	10 Sets

Description

- Disposable covers for the Hip Fix lateral supports
- Helps protect the Hip Fix lateral supports from damage or contamination
- Increases patient comfort
- Latex and phthalate free
- Single-use aids infection control

Application Area

- OR
- Anaesthesia

Measurements

- Length 205mm/8.07"
- Width 150mm/5.91"
- Height 60mm/2.36"

1x large pad cover / 2x small pad covers

Material Composition

- White foam: Polyurethane
- Blue non-woven (waterproof): 55% Woodpulp 45% Polyester
- Latex & phthalate free

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 10 sets is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- If incorrectly applied to the Hip Fix lateral supports, this product can cause pressure marks or ulcers
- Do not reuse
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624126178

Labelling

Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

Medical Device | Do Not Reuse | Latex Free | Non Sterile | UK Conformity Assessed | CE Mark

www.care-surgical.com

EU REP

Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Rev.01 23/10/2025